


PI3K-IN-30

Chemical Properties

CAS No. : 2281803-22-1
 Formula: C₂₀H₂₅F₂N₇O₃
 Molecular Weight: 449.45
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PI3K-IN-30 (compound 6d) is a potent PI3K inhibitor with IC ₅₀ values of 5.1 nM for PI3K α , 136 nM for PI3K β , 30.7 nM for PI3K γ , and 8.9 nM for PI3K δ .
Targets(IC ₅₀)	Others,PI3K

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2249 mL	11.1247 mL	22.2494 mL
5 mM	0.445 mL	2.2249 mL	4.4499 mL
10 mM	0.2225 mL	1.1125 mL	2.2249 mL
50 mM	0.0445 mL	0.2225 mL	0.445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481