


3"-O-Desmethylspinorhamnoside

Chemical Properties

CAS No. :	2220243-41-2
Formula:	C33H38O15
Molecular Weight:	674.652
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3"-O-Desmethylspinorhamnoside is a natural product for research related to life sciences. The catalog number is TN5478 and the CAS number is 2220243-41-2.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4823 mL	7.4113 mL	14.8225 mL
5 mM	0.2965 mL	1.4823 mL	2.9645 mL
10 mM	0.1482 mL	0.7411 mL	1.4823 mL
50 mM	0.0296 mL	0.1482 mL	0.2965 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481