


(E)-Guanabenz

Chemical Properties

CAS No. :	60329-03-5
Formula:	C ₈ H ₈ Cl ₂ N ₄
Molecular Weight:	231.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(E)-Guanabenz ((E)-Wy-8678), an oral central α_2 -adrenoceptor agonist, exhibits antihypertensive and neuroprotective activities by stimulating central α_2 -adrenoceptors to reduce sympathetic outflow and directly inhibiting GADD34. It is used in researching hypertension and Parkinson's disease [1] [2].
Targets(IC50)	Others,Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3275 mL	21.6375 mL	43.2751 mL
5 mM	0.8655 mL	4.3275 mL	8.655 mL
10 mM	0.4328 mL	2.1638 mL	4.3275 mL
50 mM	0.0866 mL	0.4328 mL	0.8655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481