

## Compound T130657

## Chemical Properties


CAS No. :

Formula: C17H24O11

Molecular Weight: 404.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-[3-ethenyl-5-(methoxycarbonyl)-2-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)-3,4-dihydro-2H-pyran-4-yl]acetic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.473 mL	12.3649 mL	24.7298 mL
5 mM	0.4946 mL	2.473 mL	4.946 mL
10 mM	0.2473 mL	1.2365 mL	2.473 mL
50 mM	0.0495 mL	0.2473 mL	0.4946 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481