


U 84569

Chemical Properties

CAS No. : 130736-25-3
 Formula: C₂₅H₂₃NO₄
 Molecular Weight: 401.45
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	U-84569 is a potent inhibitor of human platelet aggregation by inhibiting cyclic AMP-dependent phosphodiesterase leading to reduced platelet phospholipase C activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.491 mL	12.4549 mL	24.9097 mL
5 mM	0.4982 mL	2.491 mL	4.9819 mL
10 mM	0.2491 mL	1.2455 mL	2.491 mL
50 mM	0.0498 mL	0.2491 mL	0.4982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481