

CJM-061

## Chemical Properties


CAS No. : 1847445-87-7

Formula: C<sub>28</sub>H<sub>32</sub>N<sub>8</sub>O<sub>4</sub>

Molecular Weight: 544.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CJM-061 is a Wee1 inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8362 mL	9.1811 mL	18.3621 mL
5 mM	0.3672 mL	1.8362 mL	3.6724 mL
10 mM	0.1836 mL	0.9181 mL	1.8362 mL
50 mM	0.0367 mL	0.1836 mL	0.3672 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481