


WIN 68577

Chemical Properties

CAS No. : 151928-33-5
 Formula: C₃₆H₄₄N₈O₉
 Molecular Weight: 732.795
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	WIN 68577 is a natural product that can be used as a reference standard. The CAS number of WIN 68577 is 151928-33-5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3646 mL	6.8231 mL	13.6463 mL
5 mM	0.2729 mL	1.3646 mL	2.7293 mL
10 mM	0.1365 mL	0.6823 mL	1.3646 mL
50 mM	0.0273 mL	0.1365 mL	0.2729 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481