

4-butyl-3-methoxy-5-methylbenzoic acid

Chemical Properties


CAS No. :

Formula: C₁₃H₁₈O₃

Molecular Weight: 222.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-butyl-3-methoxy-5-methylbenzoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4988 mL	22.4942 mL	44.9883 mL
5 mM	0.8998 mL	4.4988 mL	8.9977 mL
10 mM	0.4499 mL	2.2494 mL	4.4988 mL
50 mM	0.090 mL	0.4499 mL	0.8998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481