

Compound N039-0037

Chemical Properties


CAS No. :

Formula: C₂₅H₃₅N₀O₄

Molecular Weight: 413.558

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N039-0037 is a useful organic compound for research related to life sciences and the catalog number is T131268.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.418 mL	12.0901 mL	24.1803 mL
5 mM	0.4836 mL	2.418 mL	4.8361 mL
10 mM	0.2418 mL	1.209 mL	2.418 mL
50 mM	0.0484 mL	0.2418 mL	0.4836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481