


NVP-BHS345

Chemical Properties

CAS No. : 1798292-13-3
 Formula: C₂₅H₂₀N₆O
 Molecular Weight: 420.47
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	NVP-BHS345 is a dual inhibitor of TORC1 and TORC2.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3783 mL	11.8915 mL	23.7829 mL
5 mM	0.4757 mL	2.3783 mL	4.7566 mL
10 mM	0.2378 mL	1.1891 mL	2.3783 mL
50 mM	0.0476 mL	0.2378 mL	0.4757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481