


## Ethyl 4-piperidinecarboxylate

## Chemical Properties

CAS No. :	1126-09-6
Formula:	C <sub>8</sub> H <sub>15</sub> NO <sub>2</sub>
Molecular Weight:	157.21
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ethyl 4-piperidinecarboxylate, with CAS No. 1126-09-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Ethyl 4-piperidinecarboxylate provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.3609 mL	31.8046 mL	63.6092 mL
5 mM	1.2722 mL	6.3609 mL	12.7218 mL
10 mM	0.6361 mL	3.1805 mL	6.3609 mL
50 mM	0.1272 mL	0.6361 mL	1.2722 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481