


2-Chloro-2'-C-methyladenosine

Chemical Properties

CAS No. :	205171-11-5
Formula:	C ₁₁ H ₁₄ ClN ₅ O ₄
Molecular Weight:	315.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Chloro-2'-C-methyladenosine is a 2'-C-Methyl nucleoside; Halo-nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1675 mL	15.8373 mL	31.6746 mL
5 mM	0.6335 mL	3.1675 mL	6.3349 mL
10 mM	0.3167 mL	1.5837 mL	3.1675 mL
50 mM	0.0633 mL	0.3167 mL	0.6335 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481