

p-3-Methylamino propyl phenol

Chemical Properties


CAS No. :

Formula: C₁₀H₁₅NO

Molecular Weight: 165.236

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	p-3-Methylamino propyl phenol is a useful organic compound for research related to life sciences and the catalog number is T125046.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0518 mL	30.259 mL	60.518 mL
5 mM	1.2104 mL	6.0518 mL	12.1036 mL
10 mM	0.6052 mL	3.0259 mL	6.0518 mL
50 mM	0.121 mL	0.6052 mL	1.2104 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481