


FG-7080

## Chemical Properties

CAS No. : 100332-17-0  
 Formula: C<sub>19</sub>H<sub>22</sub>FNO<sub>2</sub>  
 Molecular Weight: 315.38  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	FG-7080 is a non-tricyclic serotonin reuptake inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1708 mL	15.8539 mL	31.7078 mL
5 mM	0.6342 mL	3.1708 mL	6.3416 mL
10 mM	0.3171 mL	1.5854 mL	3.1708 mL
50 mM	0.0634 mL	0.3171 mL	0.6342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481