

## Crosemperine

## Chemical Properties


CAS No. :

Formula: C19H30NO6

Molecular Weight: 368.449

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Crosemperine is a useful organic compound for research related to life sciences and the catalog number is T125083.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7141 mL	13.5704 mL	27.1407 mL
5 mM	0.5428 mL	2.7141 mL	5.4281 mL
10 mM	0.2714 mL	1.357 mL	2.7141 mL
50 mM	0.0543 mL	0.2714 mL	0.5428 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481