


2-Phenylethanol

Chemical Properties

CAS No. :	60-12-8
Formula:	C ₈ H ₁₀ O
Molecular Weight:	122.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Phenylethanol (Phenylethanol) is an antimicrobial, antiseptic, and disinfectant. It is used also as an aromatic essence and preservative in pharmaceuticals and perfumery.
Targets(IC50)	Antibacterial, Antibiotic, Virus Protease

Solubility Information

Solubility	H ₂ O: 20 mg/mL (163.72 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.186 mL	40.9299 mL	81.8599 mL
5 mM	1.6372 mL	8.186 mL	16.372 mL
10 mM	0.8186 mL	4.093 mL	8.186 mL
50 mM	0.1637 mL	0.8186 mL	1.6372 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lingappa BT, et al. Science. 1969 Jan 10;163(3863):192-4.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481