

1-Heptacosanol

Chemical Properties

CAS No. :	2004-39-9
Formula:	C ₂₇ H ₅₆ O
Molecular Weight:	396.73
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1-Heptacosanol is a natural product for research related to life sciences. The catalog number is TN2531 and the CAS number is 2004-39-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5206 mL	12.603 mL	25.2061 mL
5 mM	0.5041 mL	2.5206 mL	5.0412 mL
10 mM	0.2521 mL	1.2603 mL	2.5206 mL
50 mM	0.0504 mL	0.2521 mL	0.5041 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481