


SKF89976A HCl

Chemical Properties

CAS No. : 85375-85-5
 Formula: C₂₂H₂₅NO₂
 Molecular Weight: 335.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SKF89976A is a potent GABA uptake inhibitor, which is selective for GAT-1 (IC ₅₀ values are 0.13, 550, 944 and 7210 μM for hGAT-1, rGAT-2, hGAT-3 and hBGT-1 respectively). SKF89976A inhibits transport current competitively (K _i = 7 μM) and transmitter-gated current non-competitively (K _i = 0.03 nM). It is able to pass the blood-brain barrier after systemic administration and is active in vivo.
Targets(IC ₅₀)	Others,GABA Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9812 mL	14.9058 mL	29.8116 mL
5 mM	0.5962 mL	2.9812 mL	5.9623 mL
10 mM	0.2981 mL	1.4906 mL	2.9812 mL
50 mM	0.0596 mL	0.2981 mL	0.5962 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481