


Sunitinib

Chemical Properties

CAS No. :	557795-19-4
Formula:	C ₂₂ H ₂₇ N ₄ O ₂
Molecular Weight:	398.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sunitinib (SU 11248) is a multi-targeted receptor tyrosine kinase (RTK) inhibitor that inhibits VEGFR2 and PDGFR β (IC ₅₀ =80/2 nM). It exhibits antitumor activity and is used for treating kidney cancer and gastrointestinal tumors.
Targets(IC ₅₀)	Apoptosis, Mitophagy, FLT, Autophagy, c-Kit, IRE1, PDGFR, VEGFR
In vitro	<p>METHODS: FLT3-ITD cell line MV4;11 and FLT3-WT cell line OC1-AML5 were treated with Sunitinib (0.001-10 μM) for 48 h, and cell viability was assayed using Alamar blue assay.</p> <p>RESULTS: Sunitinib significantly inhibited the proliferation of MV4;11 and OC1-AML5 cells in a dose-dependent manner, with IC₅₀s of 8 nM and 14 nM, respectively. [1]</p> <p>METHODS: GIST-T1 gastrointestinal mesenchymal stromal tumor cells were treated with Sunitinib (10-40 nM) for 48 h, and the expression levels of target proteins were detected by Western Blot.</p> <p>RESULTS: Sunitinib inhibited the autophosphorylation of c-KIT in a dose-dependent manner, and effectively blocked the phosphorylation of the downstream effectors of c-KIT, Akt and ERK, without affecting the phosphorylated forms of STAT3 and STAT5. [2]</p>
In vivo	<p>METHODS: To assay antitumor activity in vivo, Sunitinib (1-40 mg/kg in a citrate-buffered solution (pH 3.5)) was administered by gavage to athymic nu/nu mice harboring human leukemia tumors MV4;11 once a day until the tumors regressed.</p> <p>RESULTS: Sunitinib demonstrated dose-dependent efficacy when administered at 40 and 20 mg/kg/d and regressed established large subcutaneous tumors. [1]</p> <p>METHODS: To assay antitumor activity in vivo, Sunitinib (40-80 mg/kg in CMC) was administered by gavage to athymic nu/nu mice injected with breast cancer cells MDA-MB-435/HAL once a day for twenty-one days.</p> <p>RESULTS: Sunitinib treatment of hormonal mice resulted in a significant reduction in serum PYD levels, confirming the in vivo inhibitory effect of osteolysis. Using BLI, Sunitinib showed 64% inhibition of bone tumor growth at 40 mg/kg/day for 21 days of treatment and 89% inhibition at 80 mg/kg/day. [3]</p>
Kinase Assay	Biochemical Tyrosine Kinase Assays: IC ₅₀ values for Sunitinib against VEGFR2 (Flk-1) and PDGFR β are determined using glutathione S-transferase fusion proteins containing the complete cytoplasmic domain of the RTK. Biochemical tyrosine kinase assays to quantitate the trans-phosphorylation activity of VEGFR2 (Flk-1) and PDGFR β are performed in 96-well microtiter plates precoated (20 μ g/well in PBS; incubated overnight at 4 °C) with the peptide substrate poly-Glu,Tyr (4:1). Excess protein binding sites are blocked with the addition of 1-5% (w/v) BSA in PBS. Purified GST-fusion

Kinase Assay	proteins are produced in baculovirus-infected insect cells. GST-VEGFR2 and GST-PDGFR β are then added to the microtiter wells in 2 × concentration kinase dilution buffer consisting of 100 mM HEPES, 50 mM NaCl, 40 μ M NaVO ₄ , and 0.02% (w/v) BSA. The final enzyme concentration for GST-VEGFR2 or GST-PDGFR β is 50 ng/mL. Twenty-five μ L of diluted Sunitinib are subsequently added to each reaction well to produce a range of inhibitor concentrations appropriate for each enzyme. The kinase reaction is initiated by the addition of different concentrations of ATP in a solution of MnCl ₂ so that the final ATP concentrations spanned the K _m for the enzyme, and the final concentration of MnCl ₂ is 10 mM. The plates are incubated for 5-15 minutes at room temperature before stopping the reaction with the addition of EDTA. The plates are then washed three times with TBST. Rabbit polyclonal antiphosphotyrosine antisera are added to the wells at a 1:10,000 dilution in TBST containing 0.5% (w/v) BSA, 0.025% (w/v) nonfat dry milk, and 100 μ M NaVO ₄ and incubated for 1 hour at 37 °C. The plates are then washed three times with TBST, followed by the addition of goat antirabbit antisera conjugated with horseradish peroxidase (1:10,000 dilution in TBST). The plates are incubated for 1 hour at 37 °C and then washed three times with TBST. The amount of phosphotyrosine in each well is quantitated after the addition of 2,2'-azino-di-[3-ethylbenzthiazoline sulfonate] as substrate.
Cell Research	Cells are starved overnight in medium containing 0.1% FBS prior to addition of Sunitinib and FL (50 ng/mL; FLT3-WT cells only). Proliferation is measured after 48 hours of culture using the Alamar Blue assay or trypan blue cell viability assays. Apoptosis is measured 24 hours after Sunitinib addition by Western blotting to detect cleavage of poly (ADP-ribose) polymerase (PARP) or levels of caspase-3. (Only for Reference)

Solubility Information

Solubility	DMSO: 25 mg/mL (62.74 mM), Sonication is recommended. H ₂ O: < 1 mg/mL (insoluble or slightly soluble), Ethanol: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1.25 mg/mL (3.14 mM), Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5096 mL	12.548 mL	25.096 mL
5 mM	0.5019 mL	2.5096 mL	5.0192 mL
10 mM	0.251 mL	1.2548 mL	2.5096 mL
50 mM	0.0502 mL	0.251 mL	0.5019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- O'Farrell AM, et al. SU11248 is a novel FLT3 tyrosine kinase inhibitor with potent activity in vitro and in vivo. *Blood*. 2003 May 1;101(9):3597-605.
- Wang C, Huang M, Lin Y, et al. ENO2-derived phosphoenolpyruvate functions as an endogenous inhibitor of HDAC1 and confers resistance to antiangiogenic therapy. *Nature Metabolism*. 2023: 1-22.
- Oeller M, Jaksch-Bogensperger H, Templin M, et al. Transcription Factors STAT3 and MYC Are Key Players of Human Platelet Lysate-Induced Cell Proliferation. *International Journal of Molecular Sciences*. 2022, 23(24): 15782.
- Ikezoe T, et al. Effect of SU11248 on gastrointestinal stromal tumor-T1 cells: enhancement of growth inhibition via inhibition of 3-kinase/Akt/mammalian target of rapamycin signaling. *Cancer Sci*. 2006 Sep;97(9):945-51.
- Murray LJ, et al. SU11248 inhibits tumor growth and CSF-1R-dependent osteolysis in an experimental breast cancer bone metastasis model. *Clin Exp Metastasis*. 2003;20(8):757-66.
- Xue K H, Jiang Y F, Bai J Y, et al. Melatonin suppresses Akt/mTOR/S6K activity, induces cell apoptosis, and synergistically inhibits cell growth with sunitinib in renal carcinoma cells via reversing Warburg effect. *Redox Report*. 2023, 28(1): 2251234.
- Liu T, Yue X, Chen X, et al. Nilotinib in combination with sunitinib renders MCL-1 for degradation and activates autophagy that overcomes sunitinib resistance in renal cell carcinoma. *Cellular Oncology*. 2024: 1-18.
- Abrams TJ, et al. *Mol Cancer Ther*, 2003, 2(10), 1011-1021.
- Yee KW, et al. *Blood*, 2004, 104(13), 4202-4209.
- Dai Y, Li H, Fan S, et al. Dimethyl fumarate promotes the degradation of HNF1B and suppresses the progression of clear cell renal cell carcinoma. *Cell Death & Disease*. 2025, 16(1): 1-12.
- Machine learning-enabled virtual screening indicates the anti-tuberculosis activity of aldoxorubicin and quarfloxin with verification by molecular docking, molecular dynamics simulations, and biological evaluations

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481