


CP-74725, sodium salt

Chemical Properties

CAS No. : 111477-48-6
 Formula: C₁₂H₁₁ClFNaO₄
 Molecular Weight: 296.65
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	CP-74725, sodium salt is a bio-active chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.371 mL	16.8549 mL	33.7098 mL
5 mM	0.6742 mL	3.371 mL	6.742 mL
10 mM	0.3371 mL	1.6855 mL	3.371 mL
50 mM	0.0674 mL	0.3371 mL	0.6742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481