

## N-Boc-PEG5-alcohol

## Chemical Properties

CAS No. : 1404111-67-6


Formula: C<sub>15</sub>H<sub>31</sub>N<sub>1</sub>O<sub>7</sub>

Molecular Weight: 337.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Boc-PEG5-alcohol is a PEG-based linker for PROTACs that joins two essential ligands, facilitating the selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one ligand binds an E3 ubiquitin ligase and the other targets the specific protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9638 mL	14.8188 mL	29.6375 mL
5 mM	0.5928 mL	2.9638 mL	5.9275 mL
10 mM	0.2964 mL	1.4819 mL	2.9638 mL
50 mM	0.0593 mL	0.2964 mL	0.5928 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481