


24-Methylenecycloartanol acetate

Chemical Properties

CAS No. :	1259-94-5
Formula:	C ₃₃ H ₅₄ O ₂
Molecular Weight:	482.793
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	24-Methylenecycloartanol acetate has antiviral and antimicrobial activities.
Targets(IC50)	HSV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0713 mL	10.3565 mL	20.7129 mL
5 mM	0.4143 mL	2.0713 mL	4.1426 mL
10 mM	0.2071 mL	1.0356 mL	2.0713 mL
50 mM	0.0414 mL	0.2071 mL	0.4143 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481