

SBI-0061869

Chemical Properties


CAS No. : 853760-12-0

Formula: C₂₅H₂₈O₆

Molecular Weight: 424.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---|
| Description | SBI-0061869 is a potentiator of metabotropic glutamate receptors. It was used against neurological and psychiatric disorders. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.3558 mL | 11.7788 mL | 23.5577 mL |
| 5 mM | 0.4712 mL | 2.3558 mL | 4.7115 mL |
| 10 mM | 0.2356 mL | 1.1779 mL | 2.3558 mL |
| 50 mM | 0.0471 mL | 0.2356 mL | 0.4712 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481