


Vatalanib succinate

Chemical Properties

CAS No. :	212142-18-2
Formula:	C ₂₄ H ₂₁ ClN ₄ O ₄
Molecular Weight:	464.91
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vatalanib succinate is a VEGFR, PDGFR- β , c-Kit, c-Fms, and aromatase inhibitor.
Targets(IC50)	Others, VEGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.151 mL	10.7548 mL	21.5095 mL
5 mM	0.4302 mL	2.151 mL	4.3019 mL
10 mM	0.2151 mL	1.0755 mL	2.151 mL
50 mM	0.043 mL	0.2151 mL	0.4302 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481