

## Cbz-N-PEG15-amine

## Chemical Properties

CAS No. :

Formula: C40H74N2O17

Molecular Weight: 855.02



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Cbz-N-PEG15-amine is a PEG-based linker for PROTACs that integrates two essential ligands, crucial for forming PROTAC molecules, enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs comprise two distinct ligands connected by a linker: one binds to an E3 ubiquitin ligase, and the other targets the protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1696 mL	5.8478 mL	11.6956 mL
5 mM	0.2339 mL	1.1696 mL	2.3391 mL
10 mM	0.117 mL	0.5848 mL	1.1696 mL
50 mM	0.0234 mL	0.117 mL	0.2339 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481