


5-ClTEP

Chemical Properties

CAS No. :	900779-63-7
Formula:	C ₁₂ H ₈ ClN ₅ O ₂
Molecular Weight:	289.68
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-ClTEP is a HIV-1 integrase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4521 mL	17.2604 mL	34.5209 mL
5 mM	0.6904 mL	3.4521 mL	6.9042 mL
10 mM	0.3452 mL	1.726 mL	3.4521 mL
50 mM	0.069 mL	0.3452 mL	0.6904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481