

Compound 0407-0030

## Chemical Properties


CAS No. :

Formula: C<sub>30</sub>H<sub>43</sub>NO<sub>5</sub>

Molecular Weight: 497.676

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0407-0030 is a useful organic compound for research related to life sciences and the catalog number is T131677.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0093 mL	10.0466 mL	20.0932 mL
5 mM	0.4019 mL	2.0093 mL	4.0186 mL
10 mM	0.2009 mL	1.0047 mL	2.0093 mL
50 mM	0.0402 mL	0.2009 mL	0.4019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481