

10-Methoxygambogenic acid

Chemical Properties

CAS No. :

Formula: C₃₉H₅₀O₉

Molecular Weight: 662.82

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	10-Methoxygambogenic acid is a useful organic compound for research related to life sciences and the catalog number is T125373.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5087 mL	7.5435 mL	15.0871 mL
5 mM	0.3017 mL	1.5087 mL	3.0174 mL
10 mM	0.1509 mL	0.7544 mL	1.5087 mL
50 mM	0.0302 mL	0.1509 mL	0.3017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481