


## 3-Deoxy-11,13-dihydroisosecotanaparholide

## Chemical Properties

CAS No. :	572924-43-7
Formula:	C <sub>15</sub> H <sub>20</sub> O <sub>4</sub>
Molecular Weight:	264.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Deoxy-11,13-dihydroisosecotanaparholide is a natural product for research related to life sciences. The catalog number is TN6365 and the CAS number is 572924-43-7.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7833 mL	18.9165 mL	37.8329 mL
5 mM	0.7567 mL	3.7833 mL	7.5666 mL
10 mM	0.3783 mL	1.8916 mL	3.7833 mL
50 mM	0.0757 mL	0.3783 mL	0.7567 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481