


Mdl 28815

## Chemical Properties

CAS No. : 104905-12-6  
 Formula: C<sub>24</sub>H<sub>47</sub>NO  
 Molecular Weight: 365.64  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Mdl 28815 is a potent inhibitor of 2,3-oxidosqualene cycloartenol & lanosterol cyclases.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7349 mL	13.6747 mL	27.3493 mL
5 mM	0.547 mL	2.7349 mL	5.4699 mL
10 mM	0.2735 mL	1.3675 mL	2.7349 mL
50 mM	0.0547 mL	0.2735 mL	0.547 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481