


Apigenin 4'-O-(2'',6''-di-O-E-p-coumaroyl)glucoside

Chemical Properties

CAS No. :	71781-79-8
Formula:	C ₃₉ H ₃₂ O ₁₄
Molecular Weight:	724.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Apigenin 4'-O-(2'',6''-di-O-E-p-coumaroyl)glucoside is a natural product for research related to life sciences. The catalog number is TN6304 and the CAS number is 71781-79-8.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.380 mL	6.8998 mL	13.7996 mL
5 mM	0.276 mL	1.380 mL	2.7599 mL
10 mM	0.138 mL	0.690 mL	1.380 mL
50 mM	0.0276 mL	0.138 mL	0.276 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481