


Cyperanic acid

Chemical Properties

CAS No. :	117869-89-3
Formula:	C ₁₅ H ₂₂ O ₄
Molecular Weight:	266.337
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cyperanic acid is a useful organic compound for research related to life sciences. The catalog number is T126274 and the CAS number is 117869-89-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7546 mL	18.773 mL	37.546 mL
5 mM	0.7509 mL	3.7546 mL	7.5092 mL
10 mM	0.3755 mL	1.8773 mL	3.7546 mL
50 mM	0.0751 mL	0.3755 mL	0.7509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481