


## [2-[(4-Methylpiperazin-1-yl)methyl]phenyl]methanol

## Chemical Properties

CAS No. :	91904-36-8
Formula:	C13H20N2O
Molecular Weight:	220.31
Storage:	Store under nitrogen Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	[2-[(4-Methylpiperazin-1-yl)methyl]phenyl]methanol, with CAS No. 91904-36-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. [2-[(4-Methylpiperazin-1-yl)methyl]phenyl]methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5391 mL	22.6953 mL	45.3906 mL
5 mM	0.9078 mL	4.5391 mL	9.0781 mL
10 mM	0.4539 mL	2.2695 mL	4.5391 mL
50 mM	0.0908 mL	0.4539 mL	0.9078 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481