

H-D-Allo-Ile-OH

Chemical Properties

CAS No. : 1509-35-9


Formula: C₆H₁₃NO₂

Molecular Weight: 131.17

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	H-D-Allo-Ile-OH inhibits 3[H] gabapentin binding/voltage-dependent L-type calcium channels and can be used for studying neurological and psychiatric disorders.
Targets(IC50)	Calcium Channel,Amino Acids and Derivatives

Solubility Information

Solubility	H ₂ O: 20 mg/mL (152.47 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.6237 mL	38.1185 mL	76.2369 mL
5 mM	1.5247 mL	7.6237 mL	15.2474 mL
10 mM	0.7624 mL	3.8118 mL	7.6237 mL
50 mM	0.1525 mL	0.7624 mL	1.5247 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-1144.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481