


## Coatline B

### Chemical Properties

CAS No. : 87441-89-2  
 Formula: C<sub>21</sub>H<sub>24</sub>O<sub>11</sub>  
 Molecular Weight: 452.412  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Coatline B is a useful organic compound for research related to life sciences. The catalog number is T124030 and the CAS number is 87441-89-2.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2104 mL	11.0519 mL	22.1038 mL
5 mM	0.4421 mL	2.2104 mL	4.4208 mL
10 mM	0.221 mL	1.1052 mL	2.2104 mL
50 mM	0.0442 mL	0.221 mL	0.4421 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481