


1-(2-Methoxyphenyl)piperazine

Chemical Properties

CAS No. :	35386-24-4
Formula:	C ₁₁ H ₁₆ N ₂ O
Molecular Weight:	192.262
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-(2-Methoxyphenyl)piperazine, with CAS No. 35386-24-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-(2-Methoxyphenyl)piperazine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2013 mL	26.0064 mL	52.0129 mL
5 mM	1.0403 mL	5.2013 mL	10.4026 mL
10 mM	0.5201 mL	2.6006 mL	5.2013 mL
50 mM	0.104 mL	0.5201 mL	1.0403 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481