

meso-Hannokinol

Chemical Properties


CAS No. : 79055-11-1

Formula: C₁₉H₂₄O₄

Molecular Weight: 316.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	meso-Hannokinol is a natural product for research related to life sciences. The catalog number is TN6358 and the CAS number is 79055-11-1.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1607 mL	15.8033 mL	31.6066 mL
5 mM	0.6321 mL	3.1607 mL	6.3213 mL
10 mM	0.3161 mL	1.5803 mL	3.1607 mL
50 mM	0.0632 mL	0.3161 mL	0.6321 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481