


AG-1478

Chemical Properties

CAS No. : 153436-53-4
 Formula: C₁₆H₁₄ClN₃O₂
 Molecular Weight: 315.75
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	AG-1478 (NSC-693255) (Tyrphostin AG-1478) is a selective EGFR inhibitor.
Targets(IC50)	EGFR,HCV Protease,Influenza Virus,PDGFR
In vitro	Co-administration of 0.4 mg AG-1478 with a single dose of 25 μ Ci 90Y-CHX-A''-DTPA-hu3S193 results in significantly enhanced potency compared to the drugs administered separately. AG-1478 blocks phosphorylation of EGFR at tumor sites and inhibits growth in xenograft models of A431 cells overexpressing wt EGFR and glioma expressing de2-7 EGFR. Even subtherapeutic doses of AG-1478 significantly increase the potency of cytotoxic drugs. The combination of AG-1478 and temozolomide shows synergistic antitumor activity against human glioma xenografts.
In vivo	At a concentration of 0.25 μ M, AG-1478 effectively inhibits Ang II, Ca ²⁺ ionophore, and EGF-induced MAPK activation in VSMCs, without affecting fosinopril or platelet-derived growth factor-BB (PDGF-BB)-induced MAPK activation. AG-1478 also suppresses EGF-induced mitosis in BaF/ERX and LIM1215 cells, with IC ₅₀ values of 0.07 μ M and 0.2 μ M, respectively. Furthermore, AG-1478 inhibits the function of ABC (ATP-binding cassette) transport proteins, such as ABCB1 and ABCG2, showing a more significant effect on ABCG2. In comparison to cells expressing endogenous wild-type EGFR or overexpressing exogenous wild-type EGFR (with IC ₅₀ values of 34.6 μ M and 48.4 μ M respectively), AG-1478 preferentially inhibits U87 mg cells expressing Δ EGFR, with an IC ₅₀ of 8.7 μ M. It also preferentially inhibits tyrosine kinase activity and autophosphorylation of Δ EGFR over endogenous or overexpressed exogenous wild-type EGFR.
Cell Research	Cells are exposed to different concentrations of AG-1478 for 72 hours in 96-well plates. The effects of AG-1478 on cell growth are examined using an Alamar Blue assay. A 20- μ L aliquot of Alamar Blue is added to each well, and its absorbance is determined using a Spectromax Scanning Micro plate Reader. The effects of AG-1478 are expressed as percentage of growth inhibition using untreated cells as the control (0% inhibition). Cellular DNA synthesis is determined using a [³ H]thymidine incorporation assay. (Only for Reference)

Solubility Information

A DRUG SCREENING EXPERT

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble) H2O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 10 mg/mL (31.67 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (3.17 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1671 mL	15.8353 mL	31.6706 mL
5 mM	0.6334 mL	3.1671 mL	6.3341 mL
10 mM	0.3167 mL	1.5835 mL	3.1671 mL
50 mM	0.0633 mL	0.3167 mL	0.6334 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Levitzki A, et al. Science, 1995, 267(5205), 1782-1788.

Li S T, Young T H, Huang T W. Poly (ethylene-co-vinyl alcohol) is a suitable substrate for human olfactory neuroepithelial cell differentiation in vitro through a defined regulatory pathway. Acta Biomaterialia. 2018 Mar 1; 68:204-213

Yao G, Wang X, Wang Y, et al. miR-4792 regulates inflammatory responses in Cryptococcus neoformans-infected microglia. International Journal of Molecular Medicine. 2021, 48(5): 1-15.

Han Y, et al. Cancer Res, 1996, 56(17), 3859-3861.

Eguchi S, et al. J Biol Chem, 1998, 273(15), 8890-8896.

Johns TG, et al. Proc Natl Acad Sci U S A, 2003, 100(26), 15871-15876.

Lee FT, et al. Clin Cancer Res, 2005, 11(19 Pt 2), 7080s-7086s.

Akhtar S, et al. Cationic Polyamidoamine Dendrimers as Modulators of EGFR Signaling In Vitro and In Vivo. PLoS One. 2015 Jul 13;10(7):e0132215.

Li W, et al. EGFR Inhibition Blocks Palmitic Acid-induced inflammation in cardiomyocytes and Prevents Hyperlipidemia-induced Cardiac Injury in Mice. Sci Rep. 2016 Apr 18;6:24580.

Li S T, Young T H, Huang T W. Poly (ethylene-co-vinyl alcohol) is a suitable substrate for human olfactory neuroepithelial cell differentiation in vitro through a defined regulatory pathway[J]. Acta biomaterialia. 2018 Mar 1;68:204-213.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481