


## 3-Aminopropionitrile fumarate (2:1)

## Chemical Properties

CAS No. :	2079-89-2
Formula:	C <sub>4</sub> H <sub>4</sub> O <sub>4</sub> .2C <sub>3</sub> H <sub>6</sub> N <sub>2</sub>
Molecular Weight:	256.26
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Aminopropionitrile fumarate (2:1) [ $\beta$ -Aminopropionitrile fumarate] is an organic compound and antirheumatic agent utilized in veterinary medicine.
Targets(IC50)	MMP,Others

## Solubility Information

Solubility	DMSO: 12 mg/mL (46.83 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9023 mL	19.5114 mL	39.0229 mL
5 mM	0.7805 mL	3.9023 mL	7.8046 mL
10 mM	0.3902 mL	1.9511 mL	3.9023 mL
50 mM	0.078 mL	0.3902 mL	0.7805 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Jia LX, et al. J Pathol. 2015 Jul;236(3):373-83.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481