


Dmp 811

Chemical Properties

CAS No. : 139964-19-5
 Formula: C₂₃H₂₄N₆O₂
 Molecular Weight: 416.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dmp 811 is an angiotensin II receptor antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4011 mL	12.0054 mL	24.0108 mL
5 mM	0.4802 mL	2.4011 mL	4.8022 mL
10 mM	0.2401 mL	1.2005 mL	2.4011 mL
50 mM	0.048 mL	0.2401 mL	0.4802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481