

## Ac-IEVDIDVE (TFA)

## Chemical Properties

CAS No. :

Formula: C44H69F3N8O20

Molecular Weight: 1087.06

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Ac-IEVDIDVE TFA is a short peptide sequence.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9199 mL	4.5996 mL	9.1991 mL
5 mM	0.184 mL	0.9199 mL	1.8398 mL
10 mM	0.092 mL	0.460 mL	0.9199 mL
50 mM	0.0184 mL	0.092 mL	0.184 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481