


Divinyl glycol

Chemical Properties

CAS No. :	1069-23-4
Formula:	C ₆ H ₁₀ O ₂
Molecular Weight:	114.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Divinyl glycol is an agent of biochemical.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.7612 mL	43.8059 mL	87.6117 mL
5 mM	1.7522 mL	8.7612 mL	17.5223 mL
10 mM	0.8761 mL	4.3806 mL	8.7612 mL
50 mM	0.1752 mL	0.8761 mL	1.7522 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Topa M, Ortyl J, Chachaj-Brekiesz A, Kamińska-Borek I, Pilch M, Popielarz R. Applicability of samarium(III) complexes for the role of luminescent molecular sensors for monitoring progress of photopolymerization processes and control of the thickness of polymer coatings. *Spectrochim Acta A Mol Biomol Spectrosc.* 2018 Jun 15; 199:430-440. doi: 10.1016/j.saa.2018.03.050. Epub 2018 Mar 19. PubMed PMID: 29649679.
- Rojewska M, Bartkowiak A, Strzemiecka B, Jamrozik A, Voelkel A, Prochaska K. Surface properties and surface free energy of cellulosic etc mucoadhesive polymers. *Carbohydr Polym.* 2017 Sep 1;171:152-162. doi: 10.1016/j.carbpol.2017.05.019. Epub 2017 May 5. PubMed PMID: 28578949.
- Wang H, Cheng F, He W, Zhu J, Cheng G, Qu J. Poly(ethylene) glycol hydrogel based on oxa-Michael reaction: Precursor synthesis and hydrogel formation. *Biointerphases.* 2017 Jun 1;12(2):02C414. doi: 10.1116/1.4984305. PubMed PMID: 28571325; PubMed Central PMCID: PMC5453855.
- Choi J, Joo M, Seong H, Pak K, Park H, Park CW, Im SG. Flexible, Low-Power Thin-Film Transistors Made of Vapor-Phase Synthesized High-k, Ultrathin Polymer Gate Dielectrics. *ACS Appl Mater Interfaces.* 2017 Jun 21;9(24):20808-20817. doi: 10.1021/acsami.7b03537. Epub 2017 Jun 12. PubMed PMID: 28569054.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481