


## Valifenalate

## Chemical Properties

CAS No. :	283159-90-0
Formula:	C <sub>19</sub> H <sub>27</sub> ClN <sub>2</sub> O <sub>5</sub>
Molecular Weight:	398.88
Storage:	Store at low temperature, Store under nitrogen Powder: -20°C for 3 years   In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



## Biological Description

Description	Valifenalate (IR5885) is antifungal agent, effective against various types of mildew
Targets(IC50)	Antibacterial, Antifungal

## Solubility Information

Solubility	DMSO: 41 mg/mL (102.79 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.01 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.507 mL	12.5351 mL	25.0702 mL
5 mM	0.5014 mL	2.507 mL	5.014 mL
10 mM	0.2507 mL	1.2535 mL	2.507 mL
50 mM	0.0501 mL	0.2507 mL	0.5014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481