

CS-0045569

Chemical Properties


CAS No. : 1268024-07-2

Formula: C₁₂H₁₂O₆

Molecular Weight: 252.22

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CS-0045569, also known ZUN24972, is an important chemical intermediate for drug synthesis. It was reported in WO 2020086857 (compound 27, page 57) as an intermediate for synthesis of WDR5 inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9648 mL	19.824 mL	39.6479 mL
5 mM	0.793 mL	3.9648 mL	7.9296 mL
10 mM	0.3965 mL	1.9824 mL	3.9648 mL
50 mM	0.0793 mL	0.3965 mL	0.793 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481