

2-Methylheptanoic Acid

Chemical Properties


CAS No. : 1188-02-9

Formula: C₈H₁₆O₂

Molecular Weight: 144.21

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Methylheptanoic Acid is a flavouring ingredient.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 1.45 mg/mL (10.05 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (6.93 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9343 mL	34.6717 mL	69.3433 mL
5 mM	1.3869 mL	6.9343 mL	13.8687 mL
10 mM	0.6934 mL	3.4672 mL	6.9343 mL
50 mM	0.1387 mL	0.6934 mL	1.3869 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481