

## RGD peptide (GRGDNP) (TFA) (114681-65-1 free base)

## Chemical Properties

CAS No. :

Formula: C25H39F3N10O12

Molecular Weight: 728.63

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RGD peptide (GRGDNP) (TFA) promote apoptosis through activation of conformation changes that enhance pro-caspase-3 activation and autoprocessing.
-------------	---

## Solubility Information

Solubility	DMSO: 50 mg/mL (68.62 mM),Sonication is recommended. H2O: 50 mg/mL (68.62 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3724 mL	6.8622 mL	13.7244 mL
5 mM	0.2745 mL	1.3724 mL	2.7449 mL
10 mM	0.1372 mL	0.6862 mL	1.3724 mL
50 mM	0.0274 mL	0.1372 mL	0.2745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Ruoslahti, E., and Öbrink, B. Common principles in cell adhesion Experimental Cell Research 227, 1-11 (1996).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481