


Siponimod hemifumarate

Chemical Properties

CAS No. :	1234627-85-0
Formula:	C62H74F6N4O10
Molecular Weight:	1149.282
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Siponimod (BAF-312) hemifumarate, a selective and potent modulator of sphingosine-1-phosphate (S1P) receptors, can be used in the study of multiple sclerosis (MS).
Targets(IC50)	Others,LPL Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8701 mL	4.3505 mL	8.7011 mL
5 mM	0.174 mL	0.8701 mL	1.7402 mL
10 mM	0.087 mL	0.4351 mL	0.8701 mL
50 mM	0.0174 mL	0.087 mL	0.174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481