


SCH54292

## Chemical Properties

CAS No. : 188480-51-5  
 Formula: C<sub>24</sub>H<sub>28</sub>N<sub>2</sub>O<sub>9</sub>S  
 Molecular Weight: 520.55  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	SCH-54292 is a GDP exchange inhibitor.
Targets(IC50)	Ras

## Solubility Information

Solubility	DMSO: 55 mg/mL (105.66 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.921 mL	9.6052 mL	19.2105 mL
5 mM	0.3842 mL	1.921 mL	3.8421 mL
10 mM	0.1921 mL	0.9605 mL	1.921 mL
50 mM	0.0384 mL	0.1921 mL	0.3842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Ganguly A K , Wang Y S , Pramanik B N , et al. Interaction of a novel GDP exchange inhibitor with the Ras protein.[J]. Biochemistry, 1998, 37(45):15631-7.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481