

HAEGTFTSDVSSYLE

Chemical Properties

CAS No. :


Formula: C71H103N17O28

Molecular Weight: 1642.71

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	HAEGTFTSDVSSYLE is a peptide. GLP-I analog contains the sequence.
Targets(IC50)	Others,Glucagon Receptor

Solubility Information

Solubility	H2O: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6088 mL	3.0438 mL	6.0875 mL
5 mM	0.1218 mL	0.6088 mL	1.2175 mL
10 mM	0.0609 mL	0.3044 mL	0.6088 mL
50 mM	0.0122 mL	0.0609 mL	0.1218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

CN 102920658 B

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481