

5(6)-CR110, SE [5-(and 6)-Carboxyrhodamine 110, succinimidyl ester] \*Mixed isomers\*

## Chemical Properties

CAS No. : 254732-34-8  
 Formula: C<sub>25</sub>H<sub>18</sub>ClN<sub>3</sub>O<sub>7</sub>  
 Molecular Weight: 523.92  
 Storage: Keep away from direct sunlight  
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.

## Biological Description

Description	5(6)-CR110, SE [5-(and 6)-Carboxyrhodamine 110, succinimidyl ester] *Mixed isomers* is a Fluorescein for peptide and oligonucleotide labeling. Compared to fluorescein labeling reagents such as FITC and FAM, CR110 reagents give more photostable and pH-independent bioconjugates that have the absorption maximum around the preferred 488 nm excitation wavelength. They are photostable alternative reagents superior to FITC and FAM.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9087 mL	9.5434 mL	19.0869 mL
5 mM	0.3817 mL	1.9087 mL	3.8174 mL
10 mM	0.1909 mL	0.9543 mL	1.9087 mL
50 mM	0.0382 mL	0.1909 mL	0.3817 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481